

photobiologySensors Version 0.1.6

User Guide

Pedro J. Aphalo

June 7, 2014

1 Introduction

Here we give some examples of how one can *approximately* assess the errors involved in using a broadband sensor calibrated under a different light source, and how to *approximately* correct for this discrepancy when the spectral response of the sensor, and the spectrum of both the radiation source used for calibration and the radiation source being measured are known.

```
library(photobiology)

## Loading required package: data.table
## Loading required package: lubridate
##
## Attaching package: 'lubridate'
##
## The following objects are masked from 'package:data.table':
##
##   hour, mday, month, quarter, wday, week, yday, year
## Warning: replacing previous import by 'lubridate::hour' when loading 'photobiology'
## Warning: replacing previous import by 'lubridate::mday' when loading 'photobiology'
## Warning: replacing previous import by 'lubridate::month' when loading 'photobiology'
## Warning: replacing previous import by 'lubridate::quarter' when loading 'photobiology'
## Warning: replacing previous import by 'lubridate::wday' when loading 'photobiology'
## Warning: replacing previous import by 'lubridate::week' when loading 'photobiology'
## Warning: replacing previous import by 'lubridate::yday' when loading 'photobiology'
## Warning: replacing previous import by 'lubridate::year' when loading 'photobiology'

library(photobiologyUV)
library(photobiologyLamps)
library(photobiologyFilters)
library(photobiologySensors)
```

2 Calculating sensor response

In a way very similar to how we calculate energy irradiance with the function `energy_irradiance`, we can calculate the sensor response with the function

`sensor_response`. Both functions expect by default spectral energy irradiance as input, but will accept spectral photon irradiance as input if the parameter `unit.in` is set to “photon”.

```
sun.known.irrad <-
  with(sun.data,
    energy_irradiance(w.length, s.e.irrad, CIE()))
sun.BW20.response <-
  with(sun.data,
    sensor_response(w.length, s.e.irrad, "Vital_BW_20"))
```

The effective UV energy irradiance according to CIE erythral action spectrum for this solar spectrum is 0.0818 W m^{-2} and the BW-20 sensor response is 0.0914.

```
tl12.known.irrad <-
  with(philips.tl12.bentham.data,
    energy_irradiance(w.length, s.e.irrad, CIE()))
tl12.BW20.response <-
  with(philips.tl12.bentham.data,
    sensor_response(w.length, s.e.irrad, "Vital_BW_20"))
```

The effective UV energy irradiance according to CIE erythral action spectrum for this emission spectrum for a Philips 40W/TL12 UV-B lamp (not filtered) is 3.27 W m^{-2} and the BW-20 sensor response is 3.126.

3 Assessing errors caused by spectral mismatch

The spectral response of broad band sensors differs to smaller or greater degree from that of the theoretical response based on the BSWF of interest. Differences in spectral response between individual sensors of the same model or type are quite frequent. Consequently, all the data included in this package should be considered only as approximate. Any precise calculation will require the characterization of the spectral response of the actual individual sensor unit of interest.

We will now estimate the error involved in the use of a Vital BW-20 sensor calibrated for sunlight, to measure Q-Panel UVB-313 lamps.

```
sun.BW20.calibration.factor <-
  with(sun.data,
    sensor_response(w.length, s.e.irrad,
      "Vital_BW_20",
      reference.irrad=sun.known.irrad))
```

```
tl12.BW20.calibration.factor <-
  with(philips.tl12.bentham.data,
    sensor_response(w.length, s.e.irrad,
      "Vital_BW_20",
      reference.irrad=tl12.known.irrad))
```

The BW-20 calibration factor for solar radiation is 1.118 while that for an unfiltered Philips 40W/TL12 UV-B lamps is 0.9558. In other words measurements of TL12 lamps with the BW-20 sensor calibrated for sunlight are 117 % of the true value.

Now we check a more realistic example, a TL12 lamp filtered with cellulose diacetate. We still consider the same BW-20 sensor calibrated for sunlight.

```
acetate.s.e.irrad <-  
  with(philips.tl12.bentham.data,  
        s.e.irrad *  
        calc_filter_multipliers(w.length, "acetate.115um.new"))  
tl12.act.known.irrad <-  
  with(philips.tl12.bentham.data,  
        energy_irradiance(w.length, acetate.s.e.irrad, CIE()))  
tl12.act.BW20.calibration.factor <-  
  with(philips.tl12.bentham.data,  
        sensor_response(w.length, acetate.s.e.irrad,  
                          "Vital_BW_20",  
                          reference.irrad=tl12.act.known.irrad))
```

In this case, measurements of a diacetate-filtered TL12 lamp with the BW-20 sensor calibrated for sunlight are 103 % of the true value.